

XRF Fusion Labware suitable for Claisse machines

Dimple Crucible

CODE NUMBER	DIMENSIONS mm					ESTIMATED WEIGHT (g)	DIAGRAM
	A	B	C	D	E		
CC1	37	32	20	5	43	26	
CC2	37	32	20	5	43	30	
CC8	36	35	22.3	5	43	60	
CC21 AD6343	36.5	32	20	5.5	43	26	
CC22 AD6343	36.5	32	20	5.5	43	31	

Flat Base Crucible

CODE NUMBER	DIMENSIONS mm					ESTIMATED WEIGHT (g)	DIAGRAM
	A	B	C	D	E		
CC3	35	35	17	3.5	40	21	
CC4	35	35	21	3.5	44	26	
CC24 AD6345	36.5	32	20	5.5	43	26	
CC25 AD6345	36.5	32	20	5.5	43	30	


Flat Base Side Dimple Crucible

CODE NUMBER	DIMENSIONS mm					ESTIMATED WEIGHT (g)	DIAGRAM
	A	B	C	D	E		
CC5	37	32	20	5	43	30	
CC5-WB	37	32	23	5	43		
CC26 AD6345	36.5	32	20	5.5	43	26	
CC27 AD6345	36.5	32	20	5.5	43	30	

Flat Base NO Side Dimple Crucible

CODE NUMBER	DIMENSIONS mm					ESTIMATED WEIGHT (g)	DIAGRAM
	A	B	C	D	E		
CC5-N	37	32	20	5	43	30	
CC23 AD6344	31.6	35	15	3.9	43	24	

Mould

CODE NUMBER	DIMENSIONS mm			ESTIMATED WEIGHT (g)	DIAGRAM
	A	B	C		
CM1	30/31	6	39	18 or 30	
CM9	30/32	6	39	21 or 27	
CM2	31/33	6	41	19 or 32	
CM21	31.5/32.6	3.5	41	31	
CM3	32/34	6	41	21 or 36	
CM13	32/36	6	45	40 or 50	
CM4	32/38	6	45	22 or 36	
CM22	34.5/35.6	3.5	44	44	
CM5	35/37	6	44	25 or 48	
CM8	37/40	6	49	28 or 48	
CM6	38/41	6	49	29 or 48	
CM10	38/40	6	48	34 or 43	
CM23	39.5/40.6	3.5	49	49	
CM7	40/42	6	49	34 or 48	