

Générateur de solutions créatives



Generating creative solutions

Dishes


Evaporating dishes may be ordered also with reinforced rim and heavy base. The additional cost for reinforcement is far outweighed by longer life and ease of handling. If covers are required, please specify.

The weights indicated are for pure platinum only.

Pt95/Au5 & Pt97/Ir3 please add 10% to the pure platinum weight.

Pt90/Rh10 & Pt90/Ir10 please add 15% to the pure platinum weight.

Flat Bottom Dish


CODE NUMBER	CAPACITY cc	DIAMETER mm	BASE mm	DEPTH mm	ESTIMATED WEIGHT (g)					DIAGRAM
					STD. FORM	RIM ONLY	BASE ONLY	RIM & BASE	STANDARD LID ONLY	
GD20	20	41	28	20	6	7	8	9	6	
GD25	25	42	28	23	8	9	10	12	10	
GD35	35	49	31	25	12	13	14	16	13	
GD40	40	51	34	26	15	17	18	20	14	
GD50	50	53	46	29	17	19	21	25	15	
GD60	60	58	36	30	20	21	23	27	19	
GD75	75	65	35	30	25	27	30	33	22	
GD100	100	70	51	33	33	35	37	40	25	
GD125	125	74	48	34	42	44	45	48	28	
GD150	150	77	60	41	50	53	55	60	33	
GD175	175	84	64	40	55	59	64	68	41	
GD200	200	90	56	43	67	71	75	80	44	
GD250	250	94	65	42	80	82	83	88	51	
GD300	300	110	72	45	94	98	100	105	55	
GD400	400	115	72	48	133	137	140	148	58	
GD500	500	124	73	50	170	179	185	200	61	
GD750	750	124	73	75	270	290	310	350	61	

Générateur de solutions créatives




Generating creative solutions

General Narrow Dish

CODE NUMBER	CAPACITY cc	TOP DIAMETER mm	BOTTOM DIAMETER mm	DEPTH mm	ESTIMATED WEIGHT (g)	DIAGRAM
GND100 AD6422	100	70	36	35	40	

General Cylindrical Dish

CODE NUMBER	TOOL NUMBER	Capacity cc	Top Diameter mm (A)	Bottom Diameter mm (B)	Depth (C)	ESTIMATED WEIGHT (gm)	DIAGRAM
						STANDARD FORM	
GCD20	AD6330	20	40	32	20	13	
GCD40	AD6331	40	50	42	25	19	
GCD60	AD6332	60	60	52	30	28	
GCD80	AD6334	80	65	57	33	34	
GCD100	AD6335	100	70	64	35	45	
GCD130	AD6336	130	75	62	38	52	
GCD150	AD6337	150	80	66	40	60	
GCD200	AD6338	200	90	72	45	75	
GCD300	AD6339	300	100	82	50	115	
GCD400	AD6340	400	110	92	55	135	
GCD500	AD6341	500	120	102	60	170	
GCD600	AD6342	600	130	106	65	195	

Générateur de solutions créatives



General Cylindrical Dish (Reinforced Rim)

CODE NUMBER	TOOL NUMBER	Capacity cc	Top Diameter mm (A)	Bottom Diameter mm (B)	Depth (C)	ESTIMATED WEIGHT (gm)	DIAGRAM
						Reinforced Rim	
GCD20R	AD6330	20	40	32	20	14.5	
GCD40R	AD6331	40	50	42	25	21	
GCD60R	AD6332	60	60	52	30	30	
GCD80R	AD6334	80	65	57	33	36.5	
GCD100R	AD6335	100	70	64	35	47.5	
GCD130R	AD6336	130	75	62	38	54.5	
GCD150R	AD6337	150	80	66	40	63	
GCD200R	AD6338	200	90	72	45	78	
GCD300R	AD6339	300	100	82	50	118	
GCD400R	AD6340	400	110	92	55	139	
GCD500R	AD6341	500	120	102	60	174	
GCD600R	AD6342	600	130	106	65	199	

Evaporation Dish

CODE NUMBER	CAPACITY cc	DIAMETER mm	DEPTH mm	ESTIMATED WEIGHT (g)	DIAGRAM
GED55	55	59	22	28	
GED80	80	59	30	32	
GED100	100	67	28	40	
GED125	125	69	32	42	
GED150	150	69	40	50	

Milk Analysis Dish

CODE NUMBER	CAPACITY cc	TOP DIAMETER mm	BOTTOM DIAMETER mm	DEPTH mm	ESTIMATED WEIGHT (g)	DIAGRAM
GMD13	13	42	32	12	12	
GMD30	30	50	40	17	16	
GMD45	45	55	40	25	22	