


## Standard Form



*Standard Form  
Crucibles*


**Available Capacity:**  
2.5-700ml (larger sizes  
made to order)

**Options:**  
Lids, reinforced rims  
and/or bases.



CODE NUMBER	CAPACITY ml	DIAMETER mm	DEPTH mm	STANDARD FORM	ESTIMATED WEIGHT (g)			
					RIM ONLY	BASE ONLY	RIM & BASE	STANDARD LID ONLY
GCS2.5	2.5	12	15	3	4	4	5	2
GCS5	5	20	22	4	5	5	6	2
GCS8	8	25	24	6	7	7	8	2
GCS10	10	25	28	7	7.5	8	9	3
GCS15	15	30	31	11	12	13	14	4
GCS20	20	33	35	15	16	17	18	5
GCS25	25	36	34	19	20	21	23	5
GCS30	30	40	36	22	24	26	27	7
GCS35	35	41	38	25	26	27	32	7
GCS40	40	42	42	30	32	34	36	7
GCS50	50	45	44	38	40	42	45	8
GCS60	60	45	47	46	48	50	54	9
GCS70	70	49	53	50	52	54	63	12
GCS80	80	51	53	51	54	63	72	14
GCS90	90	54	57	52	59	71	81	15
GCS100	100	56	57	59	68	73	90	17
GCS110	110	56	60	68	71	81	95	19
GCS120	120	64	66	84	88	90	98	22
GCS150	150	65	69	100	104	107	114	27
GCS200	200	80	65	150	159	165	180	30
GCS250	250	81	73	175	179	183	190	32
GCS700	700	107	105	300	305	315	320	50

## Low Profile Form



*Low Profile Crucibles*

**Available Capacity:**  
8-125ml (larger sizes made  
to order)

**Options:**  
Lids, reinforced rims  
and/or bases.

CODE NUMBER	CAPACITY ml	DIAMETER mm	DEPTH mm	STANDARD FORM	ESTIMATED WEIGHT (g)			
					RIM ONLY	BASE ONLY	RIM & BASE	STANDARD LID ONLY
GCL8	8	25	21	6	7	7	8	2
GCL10	10	27	23	7	7.5	8	9	3
GCL15	15	30	27	11	12	13	14	4
GCL20	20	33	28	15	16	17	18	5
GCL25	25	36	30	19	20	21	23	5
GCL30	30	40	32	22	24	26	27	7
GCL35	35	42	33	25	26	27	32	7
GCL40	40	43	34	30	32	34	36	7
GCL50	50	46	37	38	40	42	45	8
GCL60	60	48	39	46	48	50	54	9
GCL70	70	52	41	50	52	54	63	12
GCL75	75	52	44	50	53	58	67	13
GCL80	80	52	47	51	54	63	72	14
GCL90	90	57	48	52	59	71	81	15
GCL100	100	57	51	59	68	73	90	17
GCL110	110	57	54	68	71	81	95	19
GCL120	120	57	57	84	88	90	98	22
GCL125	125	58	59	97	101	105	114	25

*Reinforcing of the rim gives the crucible a more robust character enabling regular handling and reducing rim deformation.*